


How to login

The screenshot shows the Biolucida web interface with the 'Connections' tab selected. A 'Biolucida Connections' dialog box is open, displaying a list of connections with 'LRZ' selected. The 'New' button is circled in red, and an arrow points to the 'LRZ' entry. Another arrow points to the 'Save' button in the 'Connection Info' section. Below the dialog, the connection details are listed:

Connection Name: LRZ
URL: vm-138-246-234-196.cloud.mwn.de
Username: Student
Password: Histokurs2021!+


Ordnerstruktur im Browser

Alle benötigten Präparate für einen Kurstag finden Sie im entsprechenden Ordner. Sortierung wie im Lernzielkatalog!




Ordnerstruktur im Browser

Wenn Sie einen Ordner öffnen, werden alle Präparate des Kurstages angezeigt. Durch Doppelklick wählen Sie ein Präparat aus.




Erstes Mikroskopieren

Zoom



Über den Regler navigieren Sie durch den Z-Stack (Ebene) = „scharfstellen“

Übersichtsbild anzeigen



Aktivierung der Makroansicht